

Product Datasheet

Rat NAT2 ELISA Kit (orb441053)

Description	Rat N-Acetyltransferase 2 (NAT2) Sandwich ELISA Kit
Reactivity	Rat
Range	0.625-40ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P50298
Sensitivity	0.218ng/mL
Entrez	116632
Expiration Date	Please enquire.