



Product Datasheet

CD15 antibody (Biotin) (orb44109)

Description	Mouse Monoclonal to CD15.
Reactivity	Human
Conjugation	Biotin
Tested Applications	FC, IP
Immunogen	Human granulocytes
Target	CD15
Preservatives	Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgM
Clonality	Monoclonal
Purity	Purified antibody is conjugated with biotin LC-NHS ester under optimum conditions and unconjugated antibody and free biotin are removed by size-exclusion chromatography.
Dilution Range	Flow cytometry: recommended dilution: 1-4 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of human neutrophil granulocytes...



Flow cytometry surface staining pattern ...