

Product Datasheet

Human NET1 ELISA Kit (orb441092)

Description	Human Neuroepithelial Cell Transforming Gene 1 (NET1) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q7Z628
Sensitivity	0.125ng/mL
Entrez	10276
Expiration Date	Please enquire.