

Product Datasheet

Mouse NEUROG3 ELISA Kit (orb441098)

Description	Mouse Neurogenin 3 (NEUROG3) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P70661
Sensitivity	0.054ng/mL
Entrez	11925
Expiration Date	Please enquire.