

Product Datasheet

Mouse NFATC1 ELISA Kit (orb441102)

Description	Mouse Nuclear Factor Of Activated T-Cells, Cytoplasmic 1 (NFATC1) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	O88942
Sensitivity	0.055ng/mL
Entrez	18018
Expiration Date	Please enquire.