

Product Datasheet

Rat NID2 ELISA Kit (orb441128)

Description	Rat Nidogen 2 (NID2) Sandwich ELISA Kit
Reactivity	Rat
Range	78.12-5000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	B5DFC9
Sensitivity	32ng/mL
Entrez	302248
Expiration Date	Please enquire.