

Product Datasheet

Human NIS ELISA Kit (orb441132)

Description	Human Sodium Iodide Symporter (NIS) Sandwich ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q92911
Sensitivity	0.24ng/mL
Entrez	6528
Expiration Date	Please enquire.