

Product Datasheet

Human NMB ELISA Kit (orb441139)

Description	Human Neuromedin B (NMB) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P08949
Sensitivity	0.056ng/mL
Entrez	4828
Expiration Date	Please enquire.