

Product Datasheet

Rat NME1 ELISA Kit (orb441141)

Description	Rat Non Metastatic Cells 1, Protein NM23A Expressed In (NME1) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q05982
Sensitivity	0.108ng/mL
Entrez	191575
Expiration Date	Please enquire.