

Product Datasheet

Rabbit NOS2 ELISA Kit (orb441160)

Description	Rabbit Nitric Oxide Synthase 2, Inducible (NOS2) Sandwich ELISA Kit
Reactivity	Rabbit
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	O19114
Sensitivity	0.54ng/mL
Entrez	100009437
Expiration Date	Please enquire.