

Product Datasheet

Porcine NOS3 ELISA Kit (orb441163)

Description	Porcine Nitric Oxide Synthase 3, Endothelial (NOS3) Sandwich ELISA Kit
Reactivity	Porcine
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q28969
Sensitivity	5.4pg/mL
Entrez	397557
Expiration Date	Please enquire.