

Product Datasheet

Human NOX1 ELISA Kit (orb441170)

Description	Human Nicotinamide Adenine Dinucleotide Phosphate Oxidase 1 (NOX1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9Y5S8
Sensitivity	0.065ng/mL
Entrez	27035
Expiration Date	Please enquire.