

Product Datasheet

Human NPC2 ELISA Kit (orb441175)

Description	Human Niemann Pick Disease Type C2 (NPC2) Sandwich ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P61916
Sensitivity	0.31ng/mL
Entrez	10577
Expiration Date	Please enquire.