

Product Datasheet

Human NPHN ELISA Kit (orb441178)

Description	Human Nephrin (NPHN) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	urine, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	urine, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	O60500
Sensitivity	0.123ng/mL
Entrez	4868
Expiration Date	Please enquire.