

Product Datasheet

Mouse NPHN ELISA Kit (orb441179)

Description	Mouse Nephrin (NPHN) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	urine, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	urine, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q9QZS7
Sensitivity	0.123ng/mL
Entrez	54631
Expiration Date	Please enquire.