

Product Datasheet

Mouse NPNT ELISA Kit (orb441183)

Description	Mouse Nephronectin (NPNT) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q91V88
Sensitivity	0.051ng/mL
Entrez	114249
Expiration Date	Please enquire.