

Product Datasheet

Human NPR2 ELISA Kit (orb441186)

Description	Human Natriuretic Peptide Receptor 2 (NPR2) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P20594
Sensitivity	29.5pg/mL
Entrez	4882
Expiration Date	Please enquire.