

Product Datasheet

Rat NPY5R ELISA Kit (orb441209)

Description	Rat Neuropeptide Y Receptor Y5 (NPY5R) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q63634
Sensitivity	0.051ng/mL
Entrez	25340
Expiration Date	Please enquire.