

Product Datasheet

Bovine NPY ELISA Kit (orb441210)

Description	Bovine Neuropeptide Y (NPY) competitive inhibition ELISA Kit
Reactivity	Bovine
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q6RUW3
Sensitivity	8.52pg/mL
Entrez	504216
Expiration Date	Please enquire.