

Product Datasheet

Rabbit NSE ELISA Kit (orb441239)

Description	Rabbit Enolase, Neuron Specific (NSE) Sandwich ELISA Kit
Reactivity	Rabbit
Range	0.625-40ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	G1SK04
Sensitivity	0.255ng/mL
Entrez	100348989
Expiration Date	Please enquire.