

Product Datasheet

Rat NSMASE ELISA Kit (orb441242)

Description	Rat Neutral Sphingomyelinase (NSMASE) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	O35049
Sensitivity	0.068ng/mL
Entrez	94338
Expiration Date	Please enquire.