

Product Datasheet

Mouse NT3 ELISA Kit (orb441245)

Description	Mouse Neurotrophin 3 (NT3) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P20181
Sensitivity	0.058ng/mL
Entrez	18205
Expiration Date	Please enquire.