

Product Datasheet

Canine NT-ProBNP ELISA Kit (orb441263)

Description	Canine N-Terminal Pro Brain Natriuretic Peptide (NT-ProBNP) Sandwich ELISA Kit
Reactivity	Canine
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sensitivity	28pg/mL
Expiration Date	Please enquire.