

Product Datasheet

Mouse NTSR1 ELISA Kit (orb441272)

Description	Mouse Neurotensin Receptor 1, High Affinity (NTSR1) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	O88319
Sensitivity	0.066ng/mL
Entrez	18216
Expiration Date	Please enquire.