

Product Datasheet

Mouse NTXI ELISA Kit (orb441275)

Description	Mouse Cross Linked N-Telopeptide Of Type I Collagen (NTXI) competitive inhibition ELISA Kit
Reactivity	Mouse
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, urine and other biological fluids.
Sample Types	serum, plasma, urine and other biological fluids.
Sensitivity	0.96ng/mL
Expiration Date	Please enquire.