

Product Datasheet

Rat NUCB2 ELISA Kit (orb441279)

Description	Rat Nucleobindin 2 (NUCB2) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9JI85
Sensitivity	0.102ng/mL
Entrez	59295
Expiration Date	Please enquire.