

Product Datasheet

Human NUP62 ELISA Kit (orb441281)

Description	Human Nucleoporin 62kDa (NUP62) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P37198
Sensitivity	0.112ng/mL
Entrez	23636
Expiration Date	Please enquire.