

Product Datasheet

Human OAS1 ELISA Kit (orb441285)

Description	Human 2',5'-Oligoadenylate Synthetase 1 (OAS1) Sandwich ELISA Kit
Reactivity	Human
Range	4.687-300ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, cell lysates and other biological fluids.
Sample Types	serum, plasma, cell lysates and other biological fluids.
Uniprot ID	P00973
Sensitivity	1.8ng/mL
Entrez	4938
Expiration Date	Please enquire.