

Product Datasheet

Canine Osteocalcin ELISA Kit (orb441294)

Description	Canine Osteocalcin (OC) Sandwich ELISA Kit
Reactivity	Canine
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	F1PDX0
Sensitivity	0.244ng/mL
Entrez	403762
Expiration Date	Please enquire.