

## Product Datasheet

### Human OCIAD2 ELISA Kit (orb441297)

|                          |                                                                                                      |
|--------------------------|------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Human Ovarian Cancer Immuno Reactive Antigen Domain Containing Protein 2 (OCIAD2) Sandwich ELISA Kit |
| <b>Reactivity</b>        | Human                                                                                                |
| <b>Range</b>             | 0.312-20ng/mL                                                                                        |
| <b>Form/Appearance</b>   | Traditional                                                                                          |
| <b>Storage</b>           | All reagents should be kept according to the instructions on the packaging upon receipt.             |
| <b>Note</b>              | For research use only                                                                                |
| <b>Application notes</b> | tissue homogenates, cell lysates and other biological fluids.                                        |
| <b>Sample Types</b>      | tissue homogenates, cell lysates and other biological fluids.                                        |
| <b>Uniprot ID</b>        | <b>Q56VL3</b>                                                                                        |
| <b>Sensitivity</b>       | 0.112ng/mL                                                                                           |
| <b>Entrez</b>            | <b>132299</b>                                                                                        |
| <b>Expiration Date</b>   | Please enquire.                                                                                      |