

Product Datasheet

Mouse OCLN ELISA Kit (orb441299)

Description	Mouse Occludin (OCLN) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q61146
Sensitivity	0.117ng/mL
Entrez	18260
Expiration Date	Please enquire.