

Product Datasheet

Rat Osteocalcin ELISA Kit (orb441303)

Description	Rat Osteocalcin (OC) Sandwich ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P04640
Sensitivity	5.7pg/mL
Entrez	25295
Expiration Date	Please enquire.