

Product Datasheet

Rat OTC4 ELISA Kit (orb441306)

Description	Rat Octamer Binding Transcription Factor 4 (OCT4) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q6MG27
Sensitivity	0.06ng/mL
Entrez	294562
Expiration Date	Please enquire.