

Product Datasheet

Rat OGDH ELISA Kit (orb441314)

Description	Rat Oxoglutarate Dehydrogenase (OGDH) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q5XI78
Sensitivity	0.053ng/mL
Entrez	360975
Expiration Date	Please enquire.