

Product Datasheet

Human OGG1 ELISA Kit (orb441315)

Description	Human Oxoguanine Glycosylase 1 (OGG1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	015527
Sensitivity	0.063ng/mL
Entrez	4968
Expiration Date	Please enquire.