

Product Datasheet

Mouse OLFM4 ELISA Kit (orb441323)

Description	Mouse Olfactomedin 4 (OLFM4) Sandwich ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	D3YYD0
Sensitivity	13.5pg/mL
Entrez	380924
Expiration Date	Please enquire.