

Product Datasheet

Human ORM2 ELISA Kit (orb441340)

Description	Human Orosomucoid 2 (ORM2) Sandwich ELISA Kit
Reactivity	Human
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P19652
Sensitivity	27.8pg/mL
Entrez	5005
Expiration Date	Please enquire.