

Product Datasheet

Human OSGIN1 ELISA Kit (orb441344)

Description	Human Oxidative Stress Induced Growth Inhibitor 1 (OSGIN1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9UJX0
Sensitivity	0.064ng/mL
Entrez	29948
Expiration Date	Please enquire.