

Product Datasheet

Human OSMR ELISA Kit (orb441349)

Description	Human Oncostatin M Receptor (OSMR) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	urine, tissue homogenates, cell lysates and other biological fluids.
Sample Types	urine, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q99650
Sensitivity	0.058ng/mL
Entrez	9180
Expiration Date	Please enquire.