

Product Datasheet

Human P2RX7 ELISA Kit (orb441369)

Description	Human Purinergic Receptor P2X, Ligand Gated Ion Channel 7 (P2RX7) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q99572
Sensitivity	0.112ng/mL
Entrez	5027
Expiration Date	Please enquire.