

Product Datasheet

Rat P4Ha1 ELISA Kit (orb441373)

Description	Rat Prolyl-4-Hydroxylase Alpha Polypeptide I (P4Ha1) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P54001
Sensitivity	0.072ng/mL
Entrez	64475
Expiration Date	Please enquire.