

Product Datasheet

Human P4Ha2 ELISA Kit (orb441374)

Description	Human Prolyl-4-Hydroxylase Alpha Polypeptide II (P4Ha2) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	O15460
Sensitivity	0.063ng/mL
Entrez	8974
Expiration Date	Please enquire.