

Product Datasheet

Human PADI2 ELISA Kit (orb441378)

Description	Human Peptidyl Arginine Deiminase Type II (PADI2) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9Y2J8
Sensitivity	0.112ng/mL
Entrez	11240
Expiration Date	Please enquire.