

Product Datasheet

Rabbit PAI-1 ELISA Kit (orb441389)

Description	Rabbit SERPINE1 ELISA Kit
Reactivity	Rabbit
Range	1.56-100pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma.
Uniprot ID	O97588
Sensitivity	0.68pg/mL
Entrez	100009588
Expiration Date	Please enquire.