

Product Datasheet

Mouse PAI2 ELISA Kit (orb441392)

Description	Mouse Plasminogen Activator Inhibitor 2 (PAI2) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma.
Uniprot ID	P12388
Sensitivity	0.134ng/mL
Entrez	18788
Expiration Date	Please enquire.