

Product Datasheet

Mouse PAOX ELISA Kit (orb441402)

Description	Mouse Polyamine Oxidase (PAOX) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q8C0L6
Sensitivity	32pg/mL
Entrez	212503
Expiration Date	Please enquire.