

Product Datasheet

Rat PAR4 ELISA Kit (orb441413)

Description	Rat Protease Activated Receptor 4 (PAR4) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q920E0
Sensitivity	0.059ng/mL
Entrez	116498
Expiration Date	Please enquire.