

Product Datasheet

Human PARK7 ELISA Kit (orb441416)

Description	Human Parkinson Disease Protein 7 (PARK7) Sandwich ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q99497
Sensitivity	0.64ng/mL
Entrez	11315
Expiration Date	Please enquire.