

Product Datasheet

Human PCPE1 ELISA Kit (orb441442)

Description	Human Procollagen C Proteinase Enhancer 1 (PCPE1) Sandwich ELISA Kit
Reactivity	Human
Range	0.625-40ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q15113
Sensitivity	0.238ng/mL
Entrez	5118
Expiration Date	Please enquire.