

Product Datasheet

Human PCSK9 ELISA Kit (orb441448)

Description	Human Proprotein Convertase Subtilisin/Kexin Type 9 (PCSK9) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q8NBP7
Sensitivity	0.125ng/mL
Entrez	255738
Expiration Date	Please enquire.